

# AYRES T-PIECE & JACKSON REECE CIRCUITS


---

Meditech APAC offers an extensive Range of Specialty Ayres T-Piece Circuits for use in Paediatric Anaesthesia, resuscitation and/or for transport.

## Features:


- » *Fresh gas provided through Washington T-Piece located close to the patient end/face mask position*
- » *1L and 0.5L open tail Rebreather bag options, with and without scavenge options*
- » *Sterile versions available*
- » *Various configurations with and without fresh gas hose. ID's of 6mm, 7mm and 10mm available and varied lengths for clinical preference and to meet pressure limitations of various OEM anaesthetic machines*

*Customisation is possible:  
Contact us for a Consultation*



# JACKSON-REES CIRCUITS

- » *Great for recovery and transport*
- » *Low resistance and minimal 'dead space'*
- » *Available with 0.5L, 1L and 2L bag*



Product	Part No.	Box Qty
Ayres T Piece Circuit, 0.5L Bag with Gas Sampling Port	22AT0315GS	20
Ayres T Piece Circuit, 1L Bag, with Gas Sampling Port	22AT0322GS	20
Ayres T Circuit, 1.6M x 10mm Fresh Gas Line with ACGO Connector, Std Elbow, 1L Open-Tail Bag	22AT1610S1L	20
Ayres T Circuit 1.6M Std O2 Tubing Fresh Gas Line with ACGO adaptor, GS Elbow and 0.5L Bag	22AT1615GS	20
Ayres T Circuit 1.6M Std O2 Tubing Fresh Gas Line with ACGO adaptor, GS Elbow and 1L Bag	22AT1615GS1L	20
Ayres T-Piece Circuit, 2.5m x 7mm O2 tubing with ACGO Connector, expandable corrugated 15mm tubing, 0.5L bag Open Tail, No Elbow	22AT250705L	20
Ayres T-Piece Circuit, 1M extendible 15mm Tubing, 2.5M 7mm O2 Tubing, 0.5L Open Tail Bag.	22AT2507X1005L	20
Ayres T-Piece Circuit, 0.5L Open Tail Bag, NO Oxygen Tubing	22ATXXXXGSE05L	20
Ayres T-Piece Circuit, 1L Open Tail Bag, No Oxygen Tubing	22ATXXXXGSE	20
Ayres T-Piece Circuit, 2L Open Tail Bag, No Oxygen Tubing	22ATXXXXGSE2L	20
Ayres T-Piece Circuit, with scavenge, 0.5L, Packed Sterile	SU22AT0315GS	20
Ayres T-Piece Circuit, with scavenge, 1L, Packed Sterile	SU22AT0322GS	20
Jackson Rees Circuit with 30cm x 22mm tubing, 2.1M O2 tubing, with 2L open tail bag	22JR0322M2L-OT	20

*Contact our team for further information and support*

**JOHN PAXMAN**  
+61 (0) 434 361 558

**MICHAEL SCHIRMER**  
+61 (0) 422 458 810

**CUSTOMER SERVICE**  
+61 (0) 8 6162 9911