



meditech
SYSTEMS_APAC

CICO KIT

CAN'T INTUBATE/
CAN'T OXYGENATE
[CICO] KIT

*A complete 'Grab and Go' solution
for hospital-wide standardisation*

SUPPORTING EVERY BREATH



CANNULA

Cricothyroidotomy




- * Boldly labeled for easy identification in low-light conditions.
- * All printing covered with clear plastic - enables repeated cleaning and long-term protection.
- * Hard-wearing non-permeable Vinyl for easy cleaning (great for Infection Control!)

Meditech CICO KIT:
2 kits combined for CICO front of neck surgical access

Part No. 88CICO-1 (below)

Includes: Rapid O2 Device / 2 x BD Insyte Cannula / 5mL Luer Slip Syringe within a Procedure Tray (Sterile Kit)



KIT 1 - CANNULA CRICOTHYROIDOTOMY		
PART NO.	DESCRIPTION	EXPIRY DATE
88CICO-1	Meditech Cricothyroidotomy Kit: Rapid-O2, 2 x 14G BD Insyte Cannulae, 5mL Syringe	01/01/2025

www.meditechsystems.com.au

- * Waterproof Zipper
- * Tamper-proof seal available



SCALPEL

Bougie



Part No. 88CICO-2 (below)

Includes: ETIP Bent Cannulated Bougie with Airway Adapter / Size 6 cuffed ETT / Size 10 Scalpel / 10mL Syringe



- * Magnetically locks when folded

- * Plastic sleeve on each side for inserted Contents Slip
- * Update your own kit! Contents slip template available for download via: www.meditechsystems.com.au

KIT 2 - SCALPEL BOUGIE		
PART NO.	DESCRIPTION	EXPIRY DATE
220212ADVB	Meditech ETIP Adult Bent Cannulated Bougie, with Airway Adaptor	01/01/2025
0501	#10 Scalpel	01/01/2025
ET60PC	Meditech Endotracheal Tube, Standard, Cuffed, 6.0	01/01/2025
88TS-10	10mL Syringe	01/01/2025

www.meditechsystems.com.au



* Heavy duty, high quality Eyelets - hang displayed contents or closed folder as required



- * Clear sleeves on both the front and rear covers for 'personalisation' of the kit.
- * Rear sleeve can be used to add site-specific protocols, algorithms, or additional clinical information if desired.



Part Number	Description	Box Qty
88CICO	FULL CICO Kit. Red outer shell, including all contents ie. Cannula Cricothyroidotomy Kit and Scalpel Bougie Kit	each
88CICO-1	Cannula Cricothyroidotomy Kit ONLY (Refill)	each
88CICO-2	Scalpel Bougie Kit ONLY (Refill)	each
88CICO-E	Red outer shell ONLY (no contents)	each

**Kit is NOT provided with saline solution. Please add.*

* *Contact our team for further information and support*

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911

meditech
SYSTEMS_APAC

www.meditechsystems.com.au