

BOUGIE INTRODUCERS & THE NEW “E-TIP”

Meditech’s new E-TIP (*Endotracheal Tube Introducer Plus - with Airway Adapter*) allows positive location of the trachea and subsequent placement of the endotracheal tube in difficult airway situations.



Solid bougie



Vented bougie includes airway adapter, for adult and paediatric patients

Features:

- » 15mm male **airway adapter** included on cannulated products*
- » Depth markings in 10cm increments: Easy reference guide for insertion
- » Low friction coating: Easy insertion and operation
- » Atraumatic Coude tip: Reduces the potential for airway trauma. Angled tip assists introduction into the vocal chords
- » Bent version: Minimise storage space, or for use in pre-pack CICO Kits

The E-TIP enhances clinical safety by including the Airway Adapter in each kit!

- *Constant reminder that if passive oxygenation is to occur, it is to be administered via circuit or resuscitator and NOT directly connected to O2 tubing.
- *Immediate confirmation of end-tidal CO2 prior to inserting endotracheal tube.

**New
Paed Sizes
Available**

Product	Part No.	Box Qty
Adult Bougie with Coude tip, solid 15fr x 70cm	220212AD	10
Adult Bougie with Coude tip, vented, includes airway adapter 14fr x 70cm	220212V14ST	10
Adult Bougie with Coude tip, vented, includes airway adapter, bent 14fr x 70cm	220212V14B	10
Paediatric Bougie with Coude tip, solid 10fr x 70cm	220211PD	10
Paediatric Bougie with Coude tip, solid 6fr x 53.5cm	220206PD	10
Paediatric Bougie with Coude tip, vented 10fr x 70cm	220211PDV	10

Contact our team for further information and support

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911