

FRAZIER SUCTION TIPS

A range of high quality, plastic/aluminium Instruments for use in surgical procedures requiring fine, delicate suctioning.

Features:

- » Available as one-piece moulded Handles: short/transparent or long/opaque handles for clinical preference
- » All supplied with obturator to facilitate shaft moulding and clearing of the instrument of obstructions
- » Long Handles are provided with control vent for precise suction control
- » Single Use, Sterile



Product	Part No.	Box Quantity
Frazier Sucker, Short, Transparent handle, 8Fr.	GC1241A08	50
Frazier Sucker, Short, Transparent handle, 10Fr.	GC1241A10	50
Frazier Sucker, Short, Transparent handle, 12Fr.	GC1241A12	50
Frazier sucker, Long, white opaque handle, with vent 8Fr.	GC1242A08	25
Frazier sucker, Long, white opaque handle, with vent 10Fr.	GC1242A10	25
Frazier sucker, Long, white opaque handle, with vent 102r.	GC1242A12	25

Contact our team for further information and support

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911