

meditech

SYSTEMS_APAC

IV INFUSION SETS + IV EXTENSION SET

Meditech offers both IV Infusion and Extension Sets in a range of configurations.

IV INFUSION SETS


Meditech's growing range of IV Infusion Sets are high-quality, gravity fed IV lines, in varied configurations for everyday clinical use.

All IV sets are STERILE, LATEX-free and DEHP-free, and provide 20 drops/ml.

Features:

IV sets range from 105cm - 295cm in length. They incorporate such features as:

- » Vented drip chambers with micron filter¹, non-return valves², roller clamp with IV tube holder³, needleless valve ports for injection⁴, positional G-clamp⁵ and rotating male luer at the distal tip⁶.
- » Clinical placement of the needleless ports facilitates user friendly 'J-loop' fixation (as required) to the hand, wrist, arm or other, with room to move for the anaesthetist or clinician, at the distal injection port.



NEW IN RANGE

IV EXTENSION SET

In addition to the popular Meditech IV Infusion Sets, an IV Extension Set is now available.

- » IV Extension Set features a needle free injection port, non-return valve, additional Y-extension with capped luer lock port and rotating male luer at the distal tip.

Configurations are highly customisable, so please contact our team with any requests for your specific clinical need.



IV EXTENSION SETS

55cm IV Extension Set with needle free injection port, non-return valve and additional Y-extension with capped luer lock port.

Code

Box Quantity

88IVEXT55

300

IV INFUSION SETS

265cm Primary IV Set for gravity infusion. Vented drip chamber, non-return valve, roller clamp with IV tube holder, 2 needleless valve injection ports.

Code

Box Quantity

88IV265

300

265cm Primary IV Set for gravity infusion. Vented drip chamber, non-return valve, roller clamp with IV tube holder, 2 needleless valve injection ports.

88IV265-25

25

265cm Primary IV Set for gravity infusion. G-Clamp, Vented drip chamber, roller clamp with IV tube holder, non-return valve, 2 needleless valve injection ports.

88IV265C

300

295cm Primary IV Set for gravity infusion. Vented drip chamber, non-return valve, roller clamp with IV tube holder, 2 needleless valve injection ports. Female Luer Y-extension with pinch clamp.

88IV295-1

300

105cm Secondary IV Set for gravity infusion. Vented drip chamber, roller clamp with IV tube holder.

88IVSS105

300

105cm Secondary IV Set for gravity infusion. Vented drip chamber, roller clamp with IV tube holder.

88IVSS105-25

25

Contact our team for further information and support

JOHN PAXMAN

+61 (0) 434 361 558

MICHAEL SCHIRMER

+61 (0) 422 458 810

CUSTOMER SERVICE

+61 (0) 8 6162 9911