



LARGE BORE SUCTION HANDLES

These large internal bore suction handles are used for the rapid evacuation of liquids and semi-solids from the upper airway in Emergency situations.

Features:

- » *Large internal diameter assists with rapid suctioning in emergency airway clearance and minimises the chances of blockages (as seen with a standard Yankauer handle)*
- » *Rigid construction and long curvature at the tip follows the natural shape of a patients airway and complements the use of airway maintenance devices, if required.*
- » *4 lateral eyes as well as distal tip opening for maximum suction*
- » *Non-vented to provide constant suctioning at high pressure*
- » *Single use, packed sterile*

OD is 9.3mm
at distal tip



OD is 13mm
at widest point

Product	Part No.	Box Qty
Catheter, Large Bore, Regular Tip, 4 Lateral Eyes, non-vented, 28Fr (OD9.3mm), Sterile.	GC1202D9N2	200

Contact our team for further information and support

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911