

ORTHOPAEDIC SUCTION SETS

This Orthopaedic Suction Set with filter keeps the operative area visible, clean, and clear. The special filter built into the handle automatically traps solids and can be exchanged in a matter of seconds when filled.

Increased filter capacity makes the suction set ideal for hip arthroplasty and other procedures that produce a large amount of debris.

FEATURES:

- » Non-block design decreases the risk of damage to tissue and postoperative infection and complications
- » Super quick Filter exchange
- » Easy to insert one-way filter design
- » Rounded atraumatic tip design
- » Single use, sterile, ready to use
- » DEHP and latex free



Product	Part No.	Box Qty
Suction Set: 15cm curved tip + 23cm curved tip + straight tip, with extra filter	GC1200C01	50

Contact our team for further information and support

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911