

THE MEDITECH SAFETY-BAG

The Meditech Safety-Bag vastly enhances current practice, creating a safer environment for your patients and staff.

Safety-Bag Features:

- » Bag stands upright for optimal visibility
- » Translucent body allows visualisation of foreign material or blockages in the airway device
- » Comes with OR without Capnography port for monitoring CO₂
- » Micron-filter on expiratory port with bacterial viral protection for OT and PACU staff. Tested efficacy of 99.9+% for airborne contaminants
- » Can be provided **individually** packaged to meet infection control guidelines and Lot number traceability OR as **loose, bulk** packed to meet environmental/waste reduction policies



Product	Part No.	Box Qty
Safety-Bag enhanced airway attachment device <i>with</i> Capnography port - Individual Pkg	22STB100	100
Safety-Bag enhanced airway attachment device <i>without</i> Capnography - Individual Pkg	22STB101	100
Safety-Bag enhanced airway attachment device <i>with</i> Capnography port - Bulk Pkg	22STB100-B	100
Safety-Bag enhanced airway attachment device <i>without</i> Capnography - Bulk Pkg	22STB101-B	100

Contact our team for further information and support

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911