

SPECIMEN TRAPS

Meditech APAC's Mucus Specimen Traps are available in 40ml and 70ml variants for the accurate collection and measurement of samples.

Features:

- » 1ml increments and 5ml numerical increment markers for accuracy of specimen collection
- » Screw on cap for securement of sample once collected
- » Packed Sterile
- » Latex-Free



Product	Part No.	Box Qty
Specimen trap, 40ml with printed graduations; packed with extra screw cap, Sterile	GC8057A04	100
Specimen trap, 70ml with printed graduations; packed with extra screw cap, Sterile	GC8057A07	100

Contact our team for further information and support

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911