


meditech

SYSTEMS_APAC

STANDARD CIRCLE BREATHING CIRCUITS

Features:

- » *Adult and paediatric two-limb circuits in standard 1.8M and 3M lengths*
- » *Customised lengths such as 4.0M and 6M circuits for use in MRI settings*
- » *Customised with various size rebreathing bags and rebreathing limbs*
- » *Sterile Circuit Kits available*
- » *Adult 1.8m and 3m validated for 7-day usage *conditions apply*



X-Tendible™ tubing is a great space saver in the clinical environment!

LOOP TAIL BAG
NOW AVAILABLE



Product	Part No.	Box Qty
Adult Circuit 1.8M with 1.8M Blue Manual Loop Tail Bag limb, Fused 22mm M/M adaptor, 2L Bag	22P182218MM2LT	10
Adult Circuit 3M with 3M Blue Manual Loop Tail Bag Limb, Fused 22mm M/M adaptor, 2L Bag	22P302230MM2LT	10
Adult Circuit 4M with 3M Blue Manual Bag Limb, 2L Bag and 4M Gas Sampling Line	22P4022302L-MRI	10
Adult Circuit 6M X-Tendible™ with 6M X-Tendible™ Blue Manual Bag Limb, 2L Bag	22X602260B2L	10
Paediatric Circuit 1.8M with 1.8M Blue Manual Bag Limb with Fused 22mm M/M adaptor, 1L Bag	22P1815CSB1L	10
Paediatric Breathing Circuit 1.8M with 1.8M Blue Manual Bag Limb, 0.5L Bag	22P181515B05L	10
Paediatric Breathing Circuit 3M with 3M Blue Manual Bag Limb, 1L bag	22P3015CSB1L	10
Adult Sterile Circuit Kit: 3M Circle Circuit with Blue Manual bag Limb 22mm M fitting, 2L bag, Filter, catheter mount and elbow, 3M gas sample line	SU22P3022CGC3M	10
Paediatric Sterile Circuit Kit: 1.8M Circle Circuit, 1L bag, 3M gas sample line	SU22PEG181522-003	10
Adult Circuit, Anaesthetic, 6M Dual limb with 6M Blue Manual Bag Limb, 2L Bag	22PS602260B2L	10

Contact our team for further information and support

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911