

TEMPERATURE MONITORING


Meditech offers quality Single Patient Use Temperature Probes and Sensors, for continuous and accurate temperature measurement.

Our general purpose probes are carefully designed for comfortable insertion into the Oesophagus or Rectum.

Our Temperature Sensors for body surface present in 3 sizes for clinical variation.

Features:

- » *LATEX FREE*
- » *Single Patient Use*
- » *Accuracy + or - 0.1°C (in the range 25-45°C)*
- » *Temperature Probes available in sizes 9 and 12 Fr*
- » *Depth markings at 10, 20, and 30cm to guide insertion*



*



OTHER SENSORS AVAILABLE ON REQUEST:

- **Foley catheters with temperature (sterile)*
- **Sterile* **Tympanic*

- » *Temperature Sensors available in sizes 9mm, 7mm and 5mm*
- » *Compatible with all 400 series Temperature monitoring devices*

Product	Code	Box Quantity
Temperature Probe 9fr Oesophageal/rectal, single patient use, 400-Series	88TP1109	25
Temperature Probe 12fr Oesophageal/rectal, single patient use, 400-Series	88TP1112	25
Temperature Sensor 9mm for Body Surface Non-sterile, single patient use, 400-Series	88TSL-1190-95	25
Temperature Sensor 7mm for Body Surface Non-sterile, single patient use, 400-Series	88TSM-1190-95	25
Temperature Sensor 5mm for Body Surface Non-sterile, single patient use, 400-Series	88TSS-1190-95	25

Contact our team for further information and support

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911