

THORACIC DRAINAGE CATHETERS

Made of high quality PVC and available in straight or right angled version to assist in surgical drainage during an indication such as pneumothorax, hemothorax, effusion and empyema.

- » *Sterile, Single Use*
- » *Clear PVC for easy visualisation of drained product*
- » *Large numerical indicators every 5cm and smaller 1cm indicators to assist with accurate depth of insertion*
- » *Coloured Radiopaque line for positional verification by X-ray*
- » *6 x Lateral smooth-sided eyes and rounded distal tip opening to facilitate strong flow and ease of insertion*
- » *Solid graduated connector provided for secure attachment between the catheter and suction tubing*
- » *Large range of sizes for all clinical uses/need*



Product	Part No.	Box Qty
Straight:		
Catheter, Thoracic Chest Drain, 16Fr, Straight.	88TCD16S	20
Catheter, Thoracic Chest Drain, 20Fr, Straight.	88TCD20S	20
Catheter, Thoracic Chest Drain, 24Fr, Straight.	88TCD24S	20
Catheter, Thoracic Chest Drain, 28Fr, Straight.	88TCD28S	20
Catheter, Thoracic Chest Drain, 32Fr, Straight.	88TCD32S	20
Catheter, Thoracic Chest Drain, 36Fr, Straight.	88TCD36S	20
Angled:		
Catheter, Thoracic Chest Drain, 20Fr, Angled.	88TCD20A	10
Catheter, Thoracic Chest Drain, 24Fr, Angled.	88TCD24A	10
Catheter, Thoracic Chest Drain, 28Fr, Angled.	88TCD28A	10
Catheter, Thoracic Chest Drain, 32Fr, Angled.	88TCD32A	10

Contact our team for further information and support

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911