


CAPNOGUARD® BITE BLOCK

It is now becoming standard to routinely monitor Capnography on patients during endoscopic procedures.


The CAPNOGUARD® Bite Block incorporates both O2 supply and CO2 monitoring, and in addition, it features a soft silicone teeth-guard to protect the patient's teeth and expensive dental work.

Features:

- » Both O2 and CO2 ports for advanced monitoring/care
- » Soft teeth guard to prevent damage to the teeth



PARTNERED CAPNOGRAPHY ITEMS:



- » Highly versatile range, including add-on features such as a Micron Filter and strap for securement

Product	Part No.	Box Quantity
CAPNOGUARD® Bite Block with O2/CO2 ports and soft silicone teeth guard	22BB1102	Box of 20
CAPNOGUARD® Bite Block with O2/CO2 ports and soft silicone teeth guard, including securement strap	22BB1102ST	Box of 20
PARTNERED CAPNOGRAPHY ITEMS		
Adult Capnography Mask, with Micron Filter (no tubing)	22142CO2	Box of 50
Micron filter 0.45um luer lock fitting for gas sampling	22B3230	Box of 100
3M TIDYLINE fused O2/CO2 tubing	2211030CO2	Box of 50
4.5M TIDYLINE fused O2/CO2 tubing	2211045CO2	Box of 50
CO2 sampling line, 3m with Microstream™ connector and 0.2um filter	DM-7100/CO-10FTST	Box of 25

Contact our team for further information and support

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911