


## CAPNOGUARD® BITE BLOCK

It is now becoming standard to routinely monitor Capnography on patients during endoscopic procedures.


The CAPNOGUARD® Bite Block incorporates both O2 supply and CO2 monitoring. In addition, it features a soft silicone teeth-guard to protect the patient's teeth and expensive dental work.

### Features:

- » Both O2 and CO2 ports for advanced monitoring/care
- » Soft teeth guard to prevent damage to the teeth



### PARTNERED CAPNOGRAPHY ITEMS:



- » Highly versatile range, including add-on features such as a Micron Filter and strap for securement

Product	Part No.	Box Quantity
CAPNOGUARD® Bite Block with O2/CO2 ports and soft silicone teeth guard	22BB1102	Box of 20
CAPNOGUARD® Bite Block with O2/CO2 ports and soft silicone teeth guard, including securement strap	22BB1102ST	Box of 20
<b>PARTNERED CAPNOGRAPHY ITEMS</b>		
Adult Capnography Mask, with Micron Filter (no tubing)	22142CO2	Box of 50
Micron filter 0.45um luer lock fitting for gas sampling	22B3230	Box of 100
3M TIDYLINE fused O2/CO2 tubing	2211030CO2	Box of 50
4.5M TIDYLINE fused O2/CO2 tubing	2211045CO2	Box of 50
CO2 sampling line, 3m with Microstream™ connector and 0.2um filter	DM-7100/CO-10FTST	Box of 25

Contact our team for further information and support

**JOHN PAXMAN**  
+61 (0) 434 361 558

**MICHAEL SCHIRMER**  
+61 (0) 422 458 810

**CUSTOMER SERVICE**  
+61 (0) 8 6162 9911