

# meditech

SYSTEMS\_APAC


## E-BOARD CIRCUIT HOLDERS

4 types:

Adult, Adult Adjustable  
Paediatric, Paediatric Adjustable

### Features:

- » Double-ended, inserts under mattress
- » 30mm opening for Co-axial style circuits
- » Graduated 22mm Adult connection tapered to 15mm Paediatric connection
- » Adjustable arm optimises circuit load, reducing torque on patient airway
- » Heavy duty 6mm perspex
- » Etched re-order details for easy reference
- » Fits ALL circuits



Product	Part No.	Quantity
E-board, Adult/Paediatric Circuit Holder, graduated 22mm - 15mm openings, 30mm for Co-axial circuits, under mattress	22ACH	ea
Adult/Paediatric Circuit Holder, adjustable, 22mm, 15mm and 30mm openings with adjustable height arm, under mattress	22ACHA	ea
E-board, Paediatric Circuit Holder, 15mm openings, under mattress	22PCH	ea
Paediatric Circuit Holder, adjustable, 15mm openings with adjustable height arm, under mattress	22PCHA	ea

Contact our team for further information and support

**JOHN PAXMAN**  
+61 (0) 434 361 558

**MICHAEL SCHIRMER**  
+61 (0) 422 458 810

**CUSTOMER SERVICE**  
+61 (0) 8 6162 9911