

meditech
SYSTEMS_APAC


LARYNGEAL MASK AIRWAYS

Introducing our high quality Disposable Silicone Laryngeal Mask range, for all patient sizes.

DISPOSABLE STANDARD LARYNGEAL MASK

Features include:


- 100% high grade silicone
- One-piece mould
- Soft silicone on mask cushion
- No epiglottis bars
- 15mm universal termination point
- Colour coded pilot balloons & mls labelled on the pilot balloon (1)



(1)

DISPOSABLE REINFORCED LARYNGEAL MASK

- Strong steel spiral cage built into the wall
- Smooth insertion - maximum patient comfort
- Translucent termination point (European Standard)
- Easy identification of colour and sizing



*Translucent termination point on both Standard and Reinforced Laryngeal Mask to aid in the visualisation of blockages

DISPOSABLE STANDARD LARYNGEAL MASK

Product Code	Description	Colour	Max Cuff Inflation	Box Quantity
22LMA50	Disposable Standard Laryngeal Mask Size 5	Yellow	40 mls	Box of 5
22LMA40	Disposable Standard Laryngeal Mask Size 4	Red	30 mls	Box of 5
22LMA30	Disposable Standard Laryngeal Mask Size 3	Orange	20 mls	Box of 5
22LMA25	Disposable Standard Laryngeal Mask Size 2.5	Green	14 mls	Box of 5
22LMA20	Disposable Standard Laryngeal Mask Size 2	Clear	10 mls	Box of 5
22LMA15	Disposable Standard Laryngeal Mask Size 1.5	Black	7 mls	Box of 5
22LMA10	Disposable Standard Laryngeal Mask Size 1	Blue	4 mls	Box of 5

DISPOSABLE REINFORCED LARYNGEAL MASK

Product Code	Description	Colour	Max Cuff Inflation	Box Quantity
22LMRE50	Disposable Reinforced Laryngeal Mask Size 5	Yellow	40 mls	Box of 5
22LMRE40	Disposable Reinforced Laryngeal Mask Size 4	Red	30 mls	Box of 5
22LMRE30	Disposable Reinforced Laryngeal Mask Size 3	Orange	20 mls	Box of 5
22LMRE25	Disposable Reinforced Laryngeal Mask 2.5	Green	14 mls	Box of 5
22LMRE20	Disposable Reinforced Laryngeal Mask Size 2	Clear	10 mls	Box of 5
22LMRE15	Disposable Reinforced Laryngeal Mask Size 1.5	Black	7 mls	Box of 5
22LMRE10	Disposable Reinforced Laryngeal Mask Size 1	Blue	4 mls	Box of 5

Contact our team for further information and support

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911