

# Microstream Compatible Capnography

---

# YOUR CAPNOGRAPHY RANGE FOR EVERY SYSTEM

Meditech Systems Asia Pacific offers a complete range of Capnography products, compatible with monitoring devices utilising Microstream Gas Sampling technology. We specialise in innovation, offering a full array of patient sampling lines to meet all of your clinical needs.


## SINGLE PORT EMERGENCY NASAL CO<sub>2</sub> CANNULA

**DM-1000-AST**

*Description:*

Single port emergency Nasal CO<sub>2</sub> Cannula, 2.1m, Adult with Microstream-compatible connector and 0.2um filter

- » When one side of nostril has been used for Bronchoscopy, the other side of the nostril can be used for CO<sub>2</sub> monitoring




## NASAL CO<sub>2</sub> SAMPLING CANNULA

**DM-4100-AST**

*Description:*

Nasal CO<sub>2</sub> sampling cannula, 2.1m, Adult with Microstream-compatible connector and 0.2um filter




## LO-FLO NASAL/ORAL CO<sub>2</sub> SAMPLING CANNULA

**DM-4601/LO-AST**

*Description:*

Lo-Flo Nasal/Oral CO<sub>2</sub> sampling cannula, 2.1m, Adult with Microstream-compatible connector and 0.2um filter

- » Works with any side-stream monitor with low flow capability 50ml/min




## DUO FLOW O<sub>2</sub> + CO<sub>2</sub> SAMPLING CANNULA

**DM-6500-10FAST**

*Description:*

Duo Flow O<sub>2</sub> + CO<sub>2</sub> sampling cannula, 3m, Adult with Microstream-compatible connector and 0.2um filter

- » Able to sampling EtCO<sub>2</sub> and delivery of O<sub>2</sub> from both nostrils simultaneously
- » Minimize dilution of waveform when O<sub>2</sub> is delivered to both nostrils




## NASAL CO<sub>2</sub> SAMPLING CANNULA W/O<sub>2</sub> DELIVERY

**DM-4000-10FAST**

*Description:*

Nasal CO<sub>2</sub> sampling cannula with O<sub>2</sub> delivery, 3m, Adult with Microstream-compatible connector and 0.2um filter




---

## NASAL/ORAL DUAL CO<sub>2</sub> SAMPLING CANNULA WITH O<sub>2</sub> DELIVERY

**DM-4601-10FTAST**

*Description:*

Nasal/Oral Dual CO<sub>2</sub> sampling cannula with O<sub>2</sub> delivery, 3m, Adult with Microstream-compatible connector and 0.2um filter




---

## CAPNOMASK O<sub>2</sub> + CO<sub>2</sub> SAMPLING MASK

**CO2-MAS-10FTAST**

*Description:*

Capnomask O<sub>2</sub> + CO<sub>2</sub> sampling mask, 3m, Adult with Microstream-compatible connector and 0.2um filter




---

## ALL-PURPOSE SAMPLING LINE WITH T-PIECE CONNECTOR (AIRWAY)

**DM-7700-10FTST**

*Description:*

All-purpose sampling line with T-Piece Connector (airway), 3m, Adult with Microstream-compatible connector and 0.2um filter




---

## ALL-PURPOSE SAMPLING LINE WITH T-PIECE CONNECTOR WITH NAFION® TUBING (AIRWAY)

**DM-7700-10FTLT**

*Description:*

All-purpose sampling line with T-Piece Connector (airway), 3m, Adult with Microstream-compatible connector, 0.2um filter and Nafion® (dryer) tubing



---

## CO<sub>2</sub> SAMPLING LINE WITH MICROSTREAM-COMPATIBLE CONNECTOR

**DM-7100/CO-10FTST**

*Description:*

CO<sub>2</sub> sampling line 3m with Microstream-compatible connector and 0.2um filter




---

## CO<sub>2</sub> SAMPLING LINE 3M WITH MICROSTREAM-COMPATIBLE CONNECTOR, 0.2UM FILTER AND NAFION (DRYER) TUBING

**DM-7100/CO-10FTLT**

*Description:*

CO<sub>2</sub> sampling line 3m with Microstream-compatible connector and 0.2um filter and Nafion (dryer) tubing



# ORDERING INFORMATION

Product	Part No.	Box Qty
Single port emergency Nasal CO <sub>2</sub> Cannula, 2.1m, Adult with Microstream-compatible connector and 0.2um filter	DM-1000-AST	25
Nasal CO <sub>2</sub> sampling cannula, 2.1m, Adult with Microstream-compatible connector and 0.2um filter	DM-4100-AST	25
Lo-Flo Nasal/Oral CO <sub>2</sub> sampling cannula, 2.1m, Adult with Microstream-compatible connector and 0.2um filter	DM-4601/LO-AST	25
Duo Flow O <sub>2</sub> + CO <sub>2</sub> sampling cannula, 3m, Adult with Microstream-compatible connector and 0.2um filter	DM-6500-10FTAST	25
Nasal CO <sub>2</sub> sampling cannula with O <sub>2</sub> delivery, 3m, Adult with Microstream-compatible connector and 0.2um filter	DM-4000-10FTAST	25
Nasal/Oral Dual CO <sub>2</sub> sampling cannula with O <sub>2</sub> delivery, 3m, Adult with Microstream-compatible connector and 0.2um filter	DM-4601-10FTAST	25
Capnomask O <sub>2</sub> + CO <sub>2</sub> sampling mask, 3m, Adult with Microstream-compatible connector and 0.2um filter	CO2-MAS-10FTAST	25
All-purpose sampling line with T-Piece Connector (airway), 3m, Adult with Microstream-compatible connector and 0.2um filter	DM-7700-10FTST	25
All-purpose sampling line with T-Piece Connector (airway), 3m, Adult with Microstream-compatible connector, 0.2um filter and Nafion (dryer) tubing	DM-7700-10FTLT	25
CO <sub>2</sub> sampling line 3m with Microstream-compatible connector and 0.2um filter	DM-7100/CO-10FTST	25

\*All sampling lines are configurable to monitor-specific connectors (Microstream, male luer lock, female luer lock). Available in both Adult and Paediatric sizes, in 2.1m or 3m lengths.

**meditech**  
SYSTEMS\_APAC

4/71 McDonald Crescent  
Bassendean WA 6934  
P. +61 (0) 8 6162 9911  
F. +61 (0) 8 6162 9816  
E. admin@meditechsystems.com.au

SUPPORTING EVERY BREATH