

# meditech

SYSTEMS\_APAC


## NEBULISERS

Nebulisers are a drug delivery device used to administer medication in the form of a mist inhaled into the lungs.

Meditech Systems APAC offers unique 'Single Patient Use' Nebulisers that allow the re-use (with cleaning) of the product on the *same patient* during a 24hr period.

### Features:

- » Operates at various angles without loss of performance
- » 85% Particles produced are within 1-5 microns size
- » Highly efficient and even delivery of drug to patient during treatment
- » Wide range of nebulisers to suit every clinical need
- » Single Patient Use, NOT Single Use!!!



### CONSIDER THIS!

A patient on a 4th-hourly nebuliser regime uses SIX single use 'competitor' nebulisers, to Meditech's ONE!

Product	Part No.	Quantity
Medi-Mist™ Nebuliser Only	227760	Box of 100
Medi-Mist™ Nebuliser, 10mL, operates at all angles	228901	Box of 100
Medi-Mist™ Nebuliser with Adult Mask, 2.1m O2 Tubing	227762	Box of 50
Medi-Mist™ Nebuliser with Paediatric Mask, 2.1m O2 Tubing	227772	Box of 50
Medi-Mist™ Nebuliser with Angled Mouthpiece, 2.1M O2 Tubing	227722	Box of 50
Medi-Mist™ Nebuliser with T-Piece Mouthpiece, 2.1M O2 Tubing	227723	Box of 50

Contact our team for further information and support

**JOHN PAXMAN**  
+61 (0) 434 361 558

**MICHAEL SCHIRMER**  
+61 (0) 422 458 810

**CUSTOMER SERVICE**  
+61 (0) 8 6162 9911