


SPECIALTY BREATHING CIRCUITS

Meditech Systems APAC brings you a range of specialised Breathing Circuits, including *Entonox Analgesic Circuits*, *Ayres T-Piece Circuits*, *Mapleson C Circuits*, *Jackson-Rees Circuits* and *Bain Circuits*, for all your clinical needs.


Ayres T-Piece Circuits

- » Compact, with no valves
- » Highly configurable, including 7mm or 10mm fresh gas line with ACGO Connector or soft O₂ Flowmeter fitting, 0.5L or 1L open tail re-breathing bags




Entonox Analgesic Circuits

- » Ergonomic design allows circuit to drape over patient's hand - no need to 'hold' the circuit
- » Expandable to 2.8M allowing patient movement in the birthing suite
- » One-way inspiratory/expiratory flow valves removed the needs for additional bacterial viral filter
- » Available with pink scavenge limb and scavenge connector




Mapleson C Circuits

- » APL for controlled pressures during ventilation
- » Great for transporting patients
- » Available with 7mm or 10mm fresh gas line with ACGO Connector or soft O₂ Flowmeter fitting




Jackson-Rees Circuits

- » Great for recovery and transport
- » Low resistance and minimal 'dead space'
- » Available with 0.5L, 1L and 2L bag



Bain Circuits

- » A co-axial modification of the basic T-piece system, developed to facilitate scavenging of waste anaesthetic gases
- » APL for controlled pressures during ventilation



**ALL of our Specialty Circuits can be customised to your specific clinical requirements (ie. length, O2 tubing diameter, APL valves, expandable, scavenging etc.)*

STERILE CIRCUIT KITS

Adult and Paediatric versions available where the breathing circuit must pass through the sterile field.

Adult kit comes complete with 3-in-1 filter, catheter mount and elbow kit and Adult /Paed both come with sterile 3M Gas sample line.

Product	Part No.	Box Quantity
Entonox Circuit, 2.8M, expandable, no Scavenge	22ENTX282200	Box of 20
Entonox Circuit, 2.8M, with mouthpiece + 2.8M pink Scavenge Limb tubing, 19mm Scavenge connector	22ENTX282228	Box of 20
Ayers T-Piece Circuit, 2.5M 7mm fresh gas line with ACGO Connector, expandable tubing, no elbow, 0.5L bag with open tail	22AT250705L	Box of 20
Ayers T-Piece Circuit, 1.6M 7mm fresh gas line with ACGO Connector, Gas Sampling elbow, 0.5L bag	22AT1615GS	Box of 20
Ayers T-Piece Circuit, 1.6M 10mm fresh gas line with ACGO Connector, plain elbow, 1L bag with open tail	22AT1610S1L	Box of 20
Mapleson C Circuit Kit with 30cm x 22mm tubing, 1.8M O2 tubing, ACGO Connector and soft O2 flowmetre fitting, 2L bag	22MC182203APL	Box of 20
Mapleson C Circuit Kit, 1.8M O2 tubing, ACGO Connector and soft O2 flow metre fitting, 2L bag	22MC1822APL	Box of 20
Jackson Rees Circuit with 30cm x 22mm tubing, 2.1M O2 tubing, with 2L open tail bag	22JR0322M2L-OT	Box of 20

STERILE CIRCUIT KITS

Adult Sterile Circuit Kit: 3M Circle Circuit with Blue Manual bag Limb 22mm M fitting, 2L bag, Filter, catheter mount and elbow, 3M gas sample line	SU22P3022CGC3M	Box of 10
Paediatric Sterile Circuit Kit: 1.8M Circle Circuit, 1L bag, 3M gas sample line	SU22PEG181522-003	Box of 10

Contact our team for further information and support

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911