

meditech  
SYSTEMS\_APAC

# RAPID-02®



*A safer way of delivering oxygen in a  
“Can’t Intubate, Can’t Oxygenate” (CICO) scenario.*

Designed by Anaesthetists.

The Rapid-02® has been designed to meet a critical need for a safer Oxygen delivery device for use when performing emergency jet oxygenation via cannula cricothyroidotomy.

#### Features:

- » *Flow and pressure release during expiration reduces the risk of lung hyperinflation and barotrauma, and raised intrathoracic pressure that may impair cardiac output*
- » *Sensory feedback of cannula occlusion*
- » *Cost point that enables Rapid-02® to be included in every theatre for immediate use in an emergency*
- » *Simple to use*
- » *Manufactured in the UK to full CE Standards*



*Large exhaust port allows for airway pressure release*

**\*Available with full  
Cannula Crico Kit!**

Product	Part No.	Box Qty
Rapid-02® Cricothyroidotomy Insufflation Device, each	2290-001-EA	ea
Rapid-02® Cricothyroidotomy Insufflation Device, box of 5	2290-001	Box of 5
Cricothyroidotomy Kit: Rapid-02®, 2x 14G Cannula, 5ml Syringe and Procedure Tray	SU2290-002-EA	ea
Cricothyroidotomy Kit: Rapid-02®, 2x 14G Cannula, 5ml Syringe and Procedure Tray	SU2290-002	Box of 5

*Contact our team for further information and support*

**JOHN PAXMAN**  
+61 (0) 434 361 558

**MICHAEL SCHIRMER**  
+61 (0) 422 458 810

**CUSTOMER SERVICE**  
+61 (0) 8 6162 9911