


SPIRALE[®] DRUG DELIVERY SYSTEM

PROUDLY DISTRIBUTED IN AUSTRALIA BY MEDITECH SYSTEMS APAC

Spirale DDS is a collapsible, volumising chamber specially designed to deliver aerosolised micro-drug particles from an MDI canister. This device is suitable for alveolar drug deposition in patients requiring intubation and mechanical ventilation and also compatible with Non-invasive ventilation systems.

- » *Can be used with both Adult and Paediatric Patients- where minimising dead space is important.*
- » *Designed to optimise alveolar drug deposition*
- » *No interruption in breathing circuit / closed system = less risk of contamination*
- » *Increase in drug efficacy*
- » *Eliminates environmental drug pollution.*
- » *No disruption in ventilation process or waiting time to carry out blood gas analyses*



***COMPATIBLE WITH NEW MDI CANISTERS WITH INTEGRATED COUNTER**



Product	Part No.	Box Qty
Spirale DDS Box of 10	AMDN1311-005-10	10
Spirale DDS Single	AMDN1311-005-10-EA	EACH

Contact our team for further information and support

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911