

meditech

SYSTEMS_APAC

TWIN-LIMB™ BREATHING CIRCUITS

The Meditech Systems APAC Twin-Limb™ is a single tube anaesthesia breathing circuit, that is divided by a septum running the length of the tube. There are many advantages to this design;

- » *Typically weighs up to 25% less than traditional two-limb circuits which reduces torque on the patient's airway*
- » *Septum separating inspiratory and expiratory lumens allows heat transfer from warm exhaled gasses to heat inspired gasses*
- » *Versatility with sizing of this single limb accommodates both adult and pediatric patients*



Single Limb with Septum and Internal Gas Sampling Line



| Product | Part No. | Box Quantity |
|---|----------------|--------------|
| Twin-Limb Anaesthesia Circuit, 1.8m, 1.8m Blue Manual Limb, 2L Bag | 22LBC18222L | 10 |
| Twin-Limb Anaesthesia Circuit, 1.8M with internal Gas Sampling Line, 1.8M Blue Manual Limb, 2L Bag. | 22LBC18222GL | 10 |
| Twin-Limb Anaesthesia Circuit, 1.8M with internal Gas Sampling Line, 1.8M x 15mm Blue Manual Limb, 1L Bag. | 22LBC18222GL1L | 10 |
| Twin-Limb Anaesthesia Circuit, 3M, 3M Expandable Blue Manual Limb, 2L Bag | 22LBC30222L | 10 |
| Twin-Limb Anaesthesia Circuit, 3M with internal Gas Sampling Line, 3M Expandable Blue Manual Limb, 2L Bag | 22LBC30222GL | 10 |
| Twin-Limb Anaesthesia Circuit, 3.5M with internal Gas Sampling Line, 3.5M Expandable Blue Manual Limb, 2L Bag | 22LBC35222GL | 10 |

Contact our team for further information and support

JOHN PAXMAN
+61 (0) 434 361 558

MICHAEL SCHIRMER
+61 (0) 422 458 810

CUSTOMER SERVICE
+61 (0) 8 6162 9911